


**VECTOR III® Duo-Tee® Spring**



Contact:  
**GENE SHANNON**  
 P. O. Box 633  
 101 Industrial Dr.  
 Burlington, Wi 53105  
 262-767-1440  
 1-800-517-9424

**VECTOR III® CERTIFIED / LICENSED FOR:  
 Seminars/Appliance Fabrication**

- I. DUO-TEE® SPRINGS with VECTORLOY™ wire**  
U.S. Pat. # 6604943
  - II. Method of preadjusted fabrication**  
U.S. Pat. #6935859
- VECTORLOY™ WIRE with its' characteristic color identify the genuine**

Dr. Name \_\_\_\_\_

Address \_\_\_\_\_

Phone \_\_\_\_\_

City, State, Zip \_\_\_\_\_

Patients Name\* \_\_\_\_\_

\*Patient's name must be correct for inlays. Construct and deliver to the undersigned only.

LICENSE NUMBER \_\_\_\_\_ DATE \_\_\_\_\_

Signature \_\_\_\_\_ sent \_\_\_\_\_ need \_\_\_\_\_  
(date) (date)

X if needed

**Our lab has an Orthodontic consultant available as needed for pretreatment evaluation and treatment progress assistance**

**STANDARD COLORS:**


\_\_\_\_ PINK  
 \_\_\_\_ CLEAR  
 SPECIFY OTHER COLOR  
 (AT MINIMAL CHARGE) \_\_\_\_\_  
 \_\_\_\_ NAME INLAID?

**VECTOR III® Duo-Tee® Springs retainer**  
U.S. Pat. # 6604943

- Standard pre-activated retainer  
\*\*see below
- Standard + bite plane upper buccal/ segments with anterior x-bite correction- no extra charge (\$160.00)
- RX special instructions (see below)

**Special Instructions:**

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_



USE BACK FOR ADDITIONAL INSTRUCTIONS OR COMMENTS


**Duo-Tee II@**  
 \$165.00  
 Spring →



**VECTOR III® Duo-Tee® Springs** standard cuspid to cuspid anterior appliance made pre-activated is \$160.00. Any additional springs are \$5.00 each. Custom appliances can be priced pre-fabrication.

**AFTER FILLING OUT PLEASE MAKE COPY FOR YOUR RECORDS, ALLOW 2-4 WEEKS OR CALL AHEAD FOR SPECIAL CASES OR FASTER FABRICATION**

**VECTOR III® Duo-Tee® Spring**



Contact:  
**GENE SHANNON**  
 P. O. Box 633  
 101 Industrial Dr.  
 Burlington, Wi 53105  
 262-767-1440  
 1-800-517-9424

**VECTOR III® CERTIFIED / LICENSED FOR:  
 Seminars/Appliance Fabrication**

- I. DUO-TEE® SPRINGS with VECTORLOY™ wire**  
U.S. Pat. # 6604943
  - II. Method of preadjusted fabrication**  
U.S. Pat. #6935859
- VECTORLOY™ WIRE with its' characteristic color identify the genuine**

Dr. Name \_\_\_\_\_

Address \_\_\_\_\_

Phone \_\_\_\_\_

City, State, Zip \_\_\_\_\_

Patients Name\* \_\_\_\_\_

\*Patient's name must be correct for inlays. Construct and deliver to the undersigned only.

LICENSE NUMBER \_\_\_\_\_ DATE \_\_\_\_\_

Signature \_\_\_\_\_ sent \_\_\_\_\_ need \_\_\_\_\_  
(date) (date)

X if needed

**Our lab has an Orthodontic consultant available as needed for pretreatment evaluation and treatment progress assistance**

**STANDARD COLORS:**


\_\_\_\_ PINK  
 \_\_\_\_ CLEAR  
 SPECIFY OTHER COLOR  
 (AT MINIMAL CHARGE) \_\_\_\_\_  
 \_\_\_\_ NAME INLAID?

**VECTOR III® Duo-Tee® Springs retainer**  
U.S. Pat. # 6604943

- Standard pre-activated retainer  
\*\*see below
- Standard + bite plane upper buccal/ segments with anterior x-bite correction - no extra charge (\$160.00)
- RX special instructions (see below)

**Special Instructions:**

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_



USE BACK FOR ADDITIONAL INSTRUCTIONS OR COMMENTS

**Duo-Tee II@**  
 \$165.00  
 Spring →



**VECTOR III® Duo-Tee® Springs** standard cuspid to cuspid anterior appliance made pre-activated is \$160.00. Any additional springs are \$5.00 each. Custom appliances can be priced pre-fabrication.

**AFTER FILLING OUT PLEASE MAKE COPY FOR YOUR RECORDS, ALLOW 2-4 WEEKS OR CALL AHEAD FOR SPECIAL CASES OR FASTER FABRICATION**